

Item # 26800B-461, Auxiliary Lens 0.5X



Auxiliary Lens 0.5X

Using our 0.5x auxiliary lens will decrease magnification and extend working distance.

· [SPECIFICATIONS](#) · [NOTE](#)

SPECIFICATIONS

Lens	0.5x
For Bodies	DSZ Series NSW Series SPZ Series
Weight	0.15 lb

NOTE

Please note 0.5x lens cannot be used with SPZ series bodies with stand P and PF. 0.5x lens cannot be used for stands T and TL and with NSW series with stand P
